

Phototherapy Breakthrough • Phototherapy Breakthrough

DERMA-WAND™

GERMICIDAL UVC LAMP

From NATIONAL BIOLOGICAL



Creates A Bacteria-Free Environment
For The Treatment Of Skin Disorders

- ◆ Control of time, distance and angle for repeatable results
- ◆ No gels, electrodes or contact with affected area
- ◆ Portable
- ◆ Cost-effective

THE **POWER** OF LIGHT

It Really Works Wonders!

Phototherapy Breakthrough • Phototherapy Breakthrough

DERMA-WAND™

Features

- ◆ Electronic digital timer
- ◆ Portable and compact
- ◆ Distance guard to assure consistent, repeatable treatment
- ◆ UV blocking goggles and operating manual included
- ◆ Warranty:
Five (5) years, parts and labor (excluding lamps)
- ◆ Unlimited customer support



NRTL/C
Complies with
Standards
CAN/CSA STD.
C22.2 NO.125
M1984, UL 544.



DERMA-WAND™ Specifications:

- ◆ Two UVC lamps
- ◆ High UV output ballast design
- ◆ Full reflector maximizes UV output
- ◆ Control box weight —7.75 lbs.
Wand weight — 1.75 lbs.
- ◆ Electrical: 110-120V plug-in —
60 Hz AC .5 Amps with hospital-grade grounding plug

CAUTION: Federal law restricts this device to sale by or on the order of a physician.



**NATIONAL
BIOLOGICAL
CORPORATION**

1532 Enterprise Parkway
Twinsburg, OH 44087
(800) 338-5045
Fax: (330) 425-9614
E-mail: nbc@natbiocorp.com
Web: www.natbiocorp.com

Distributed By