

GeniStim 330™

Portable Electrical Stimulator



- ✓ Reduce pain
- ✓ Heal wounds
- ✓ Increase blood circulation

Wound Healing WITH RESULTS!

High voltage pulsed current (HVPC) for wound healing, increased blood flow and pain reduction.

EASY-TO-USE!

The GeniStim 330 is an easy-to-use high voltage pulsed current (HVPC) stimulator suited to a range of physiological applications. It supplies a galvanic stimulation current for wound healing, increased blood circulation and pain reduction. A convenient timer permits automatic shut-off at the end of the therapy session.

A pulsed electric current is effective in stimulating tissue to produce healing during all phases of the healing process.

It has two active leads for treating two locations at the same time. Controls are housed in a transparent cover to prevent inadvertent altering of the settings during treatment.

Specifications:

- **Output Channel:** Two active outputs, one dispersive return
- **Cycle:** Continuous or alternating
- **On Time:** 7 seconds fixed in alternating mode, or continuous
- **Wave Form:** Pulsed galvanic monophasic exponential spikes, delivered in pairs with a fixed interpulse interval of 100 microseconds
- **Treatment Timer:** 60 minutes, 15 on / 45 off, or continuous
- **Frequency:** 1 - 100 Hz, selectable, with digital display
- **Output Voltage:** 0 - 330 V, adjustable
- **Power Source:** Two 9-volt batteries or AC adapter, 120 VAC
- **Polarity:** Positive or negative
- **Time Display:** Records total usage time in minutes

Included with the GeniStim:

- HVPC stimulator
- Lead wire
- 3 alligator clips
- 2 splitter cables
- Electrode kit with two active and one dispersive electrodes
- AC adapter
- Instruction sheet
- Carry case

www.biomation.com/wound

™GeniStim 330 is a trademark of BIOMATION

BIOMATION

Authorized Distributor

335 Perth Street
P.O. Box 156
Almonte, ON K0A 1A0
Canada

Tel: 613-256-2821
Toll Free: 1-888-667-2324
Fax: 613-256-5872
E-mail: info@biomation.com/physiotherapy